

PRODUCT INFORMATION

Pour plate method Procedure

Pour plate method Procedure:

Serial Dilution: The sample is first serially diluted to ensure that the number of colonies on the plate will be countable (ideally 30–300 CFU per plate).

Inoculation: A measured volume (usually 1 mL) of the diluted sample is added to a sterile Petri dish.

Addition of Molten Agar: About 15–18 mL of molten agar (cooled to 40–45°C) is poured into the dish over the sample.

Mixing: The plate is gently swirled to mix the sample evenly with the agar.

Solidification and Incubation: The agar is allowed to solidify, and the plate is inverted and incubated at an appropriate temperature and time.

Colony Counting: After incubation, colonies are counted. Colonies grow both on the surface and within the medium; each colony represents a viable microorganism from the original sample.