

GIARDIA/CRYPTOSPORIDIUM CHEK[®], GIARDIA II and CRYPTOSPORIDIUM II

ASSAY PROTOCOL


- 1 **Sample Preparation**
For fresh or frozen specimens, add 100 μ L of specimen to 400 μ L

Diluent. For preserved specimens, no dilution or further processing is necessary.



- 2 **Add *Diluent* and fecal samples**
Add 100 μ L of *Diluent* to the negative control well and all test wells.

Transfer 50 μ L of diluted or preserved fecal specimen to test wells containing *Diluent*.



- 3 **Add *Positive Control* and incubate**

Add 1 drop (50 μ L) of *Positive Control* to positive control wells. Gently tap wells to mix. Seal and incubate for 60 minutes at room temperature.



- 4 **Wash Assay Wells**

Wash wells with diluted *Wash Solution* the specified number of times. Slap inverted plate hard on paper towels between each wash.



- 5 **Add *Conjugate***
Add 1 drop (50 μ L) of *Conjugate* to each well. Tap to mix.

Seal and incubate for 30 minutes at room temperature.



- 6 **Wash Assay Wells**
Repeat step 4.

- 7 **Add *Substrate* reagent**
Add 2 drops of *Substrate* to each well. Incubate wells at room temperature for 10 minutes. Tap wells gently to mix during incubation.



- 8 **Add *Stop Solution***
Add 1 drop of *Stop Solution*. Tap to mix. Wait 2 minutes before reading. Refer to the package insert for interpretation of results.



Distributed by:



Alere North America, LLC.
30 South Keller Road
Orlando, FL 32810 USA

Developed and Manufactured by:




2001 Kraft Drive
Blacksburg, VA 24060-6358
Made in the USA

This chart does not contain the complete instructions for use of the *GIARDIA/CRYPTOSPORIDIUM CHEK[®]*, *GIARDIA II* test, or *CRYPTOSPORIDIUM II* test. For proper use of the assay, please read the appropriate package insert.

Visual Interpretation for *GIARDIA*/CRYPTOSPORIDIUM CHEK[®], *GIARDIA* II, and *CRYPTOSPORIDIUM* II

1. Compare the color in the well with the color in the chart below.
2. Read the control wells first. The negative control well should be colorless or have a faint yellow color (< 0.150 at OD_{450} or < 0.090 at $OD_{450/620}$). The positive control well should be \geq the 2+ color (≥ 0.500 at OD_{450} or $OD_{450/620}$).
3. **Positive** results will have a color intensity \geq the 1+, 2+, or 3+ color (≥ 0.150 at OD_{450} or ≥ 0.090 at $OD_{450/620}$).
4. **Negative** results will be colorless or resemble the negative control well in intensity of color.



Refer to the appropriate Package Insert for complete details.