
● CHROMagar™
VRE

For detection of Van A / Van B
VRE. faecalis & *VRE. faecium*



Plate Reading

- VRE. *faecalis* / VRE. *faecium*
→ pink to mauve
- *E. gallinarum* / *E. casseliflavus*
→ blue or inhibited
- other bacteria
→ inhibited

Quality Control Strains

VRE. <i>faecalis</i> ATCC® 51575	mauve, small
<i>E. faecalis</i> ATCC® 29212	inhibited
<i>E. coli</i> ATCC® 25922	inhibited
<i>C. tropicalis</i> ATCC® 1369	inhibited
<i>S. aureus</i> ATCC® 25923	inhibited
<i>E. casseliflavus</i> ATCC® 700327	inhibited
<i>E. gallinarum</i> ATCC® 49573	inhibited

ATCC® is a registered trademark of the American Type Culture Collection

For detection of Van A / Van B VRE. *faecalis* & VRE. *faecium*

Background

There are two types of vancomycin resistance in *enterococci*. The first type is intrinsic resistance (mostly VanC type but also VanD, VanE, VanF etc) found in *E. gallinarum* and *E. casseliflavus*/*E. flavescentis* and demonstrates a low-level resistance to vancomycin. The second type of vancomycin resistance in *enterococci* is acquired resistance (VanA & VanB types), mostly seen in *E. faecium* and *E. faecalis*. Therefore, to avoid the spread of this resistance to more virulent pathogens (*S. aureus*, for instance) it is crucial to promptly detect the presence of any of these two species in the patient, and accurately differentiate them from other *Enterococci*.

“Knowledge of the type of resistance is critical for infection control purposes. VanA and VanB genes are transferable and can spread from organism to organism. In contrast, VanC genes are not transferable, have been associated less commonly with serious infections, and have not been associated with outbreaks” – from CDC guidelines

Vancomycin-resistant *Enterococcus* (VRE) infections are especially aggressive and have been associated with mortality rates approaching 60% to 70%.

Medium Performance

1 SIMPLE, FAST AND RELIABLE TOOL
for the direct detection of VRE strains with transmissible resistance: this is a precious help in the implementation of the appropriate control measures to prevent the spread of VRE.

2 INTENSE COLONY COLOURS
In the new CHROMagar™ VRE media, VRE. *faecalis* and VRE. *faecium* strains are easily distinguishable by the colony colour.

In the contrary, in the Classical agar for the detection of VRE (Bile Esculine Agar supplemented with vancomycin) : (I) there is no differentiation between *E. faecalis*/*E. faecium* and the other *enterococci*; (II) it often leads to false positives of other esculine hydrolysing bacteria (like *Lactococcus*, *Pediococcus*...); (III) the black “cloud” makes plate reading difficult as well as the choice of the proper colony for further confirmatory tests.

3 FLEXIBILITY
CHROMagar™ VRE is supplied with a shelf-life of about 2 years. This allows for flexibility of use, whether in an epidemic situation with many patients to screen, or whether for random surveillance of cultures.

Medium Description

Powder Base	Total	67.3 g/L
	Agar	15.0
+ Supplement (Included in the pack)	Peptones & Yeast extract	20.0
	Salts	5.0
	Chromogenic mix	27.3
	Storage at 15/30°C - pH: 6.9 ± 0.2	
	Shelf Life	2 years
	Powder form	60 mg/L
Storage at 2/8°C	Shelf Life	2 years

Usual Samples	stools
Procedure	Direct Streaking. Incubation at 37°C, 24h. Aerobic conditions.
Scientific Publications on this product: available on www.CHROMagar.com Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com	


Order References



Please use these product references when contacting your local distributor:

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