

PRODUCT INFORMATION

Product Type: Divided Petri Dishes 90mm (DD)

Cat No. DD017 - SS AGAR / XLD AGAR

Intended Use:

SS AGAR: SS Agar is recommended for isolation and differentiation of pathogenic enteric bacilli, especially *Salmonella* and *Shigella* from clinical specimens, including feces, and food samples.

XLD AGAR: *Salmonella* isolation including detection of *Salmonella* Variants.

Principles and uses:

SS AGAR

Differentiation of gram-negative bacilli on SS Agar is based on the fermentation of lactose and subsequent absorption of neutral red; gram negative bacilli which ferment lactose produce pink to red colonies. Lactose-non fermenters, such as *Salmonella* and *Shigella*, form transparent, colorless colonies. Bile salts, brilliant green, and citrates are selective agents which inhibit gram-positive bacteria and coliforms. Sodium thiosulfate, a sulfur source, and ferric ammonium citrate, an indicator, are added to enable organisms which produce H₂S to form black-centered colonies, including some strains of *Salmonella*. The medium contains a bile salts mixture to inhibit the growth of gram-positive organisms.

Limitations SS AGAR:

1. Organisms other than *Salmonella* and *Shigella* which grow on this medium may be differentiated by their ability to ferment lactose and form pink or red colonies.
2. SS Agar is not recommended for primary isolation of *Shigella* spp. because it is highly selective and may inhibit some strains.

XLD AGAR:

A selective medium for the isolation of *salmonellae* and *shigella* from clinical specimens and foods. Xylose-Lysine-Desoxycholate Agar was originally formulated by Taylor for the isolation and identification of *shigella* from stool specimens. It has since been found to be a satisfactory medium for the isolation and presumptive identification of both *salmonellae* and *shigella*. It relies on xylose fermentation, lysine decarboxylation and production of hydrogen sulphide for the primary differentiation of *shigella* and *salmonellae* from non-pathogenic bacteria.

Rapid xylose fermentation is almost universal amongst enteric bacteria, except for members of the *Shigella*, *Providencia* and *Edwardsiella* genera. Xylose is thus included in the medium so that *Shigella* spp. may be identified by a negative reaction. *Salmonella* spp. are differentiated from non-pathogenic xylose fermenters by the incorporation of lysine in the medium. *Salmonellae* exhaust the xylose and decarboxylate the lysine, thus altering the pH to alkaline and mimicking the *Shigella* reaction. However, the presence of *Salmonella* and *Edwardsiella* spp. is differentiated from that of *shigella* by a hydrogen sulphide indicator.

The high acid level produced by fermentation of lactose and sucrose, prevents lysine-positive coliforms from reverting the pH to an alkaline value, and non-pathogenic hydrogen sulphide producers do not decarboxylate lysine. The acid level also prevents blackening by these micro-organisms until after the 18–24-hour examination for pathogens.

Sodium desoxycholate is incorporated as an inhibitor in the medium. The concentration used allows for the inhibition of coliforms without decreasing the ability to support *shigella* and *salmonellae*. The recovery of *Shigella* spp. has previously been neglected despite the high incidence of shigellosis. This has been due to inadequate isolation media.

The sensitivity and selectivity of X.L.D. Agar exceeds that of the traditional plating media e.g. Eosin Methylene Blue, *Salmonella-Shigella* and Bismuth Sulphite agars, which tend to suppress the growth of *shigella*. Many favorable comparisons between X.L.D. Agar and these other media have been recorded in the literature.

The recovery of salmonellae and shigellae is not obscured by profuse growth of other species therefore X.L.D. Agar is ideal for the screening of samples containing mixed flora and suspected of harboring enteric pathogens e.g. medical specimens or food products. Chadwick, Delisle and Byer recommended the use of this medium as a diagnostic aid in the identification of *Enterobacteriaceae*.

Reference XLD AGAR:

1. Taylor W.I. (1965) Am. J. Clin. Path. 44. 471-475.
2. McCarthy M.D. (1966) N.Z. J. Med. Lab. Technol. 20. 127-131.
3. Isenberg H.D., Kominos S. and Sigel M. (1969) Appl. Microbiol . 18. 656-659.
4. Taylor W. I. and Harris B. (1965) Am. J. Clin. Path. 44. 476-479.
5. Taylor W. I. and Harris B. (1967) Am. J. Clin. Path. 48. 350-355.
6. Taylor W. I. and Schelhart D. (1967) Am. J. Clin. Path. 48. 356-362.
7. Taylor W.I. and Schelhart D. (1966) Appl. Microbiol . 16. 1387-1392.
8. Rollender M.A., Beckford O., Belsky R.D. and Kostroff B. (1969) Am. J. Clin. Path. 51. 284-286.
9. Taylor W. I. and Schelhart D. (1969) Appl. Microbiol . 18. 393-395.
10. Dunn C. and Martin W.J. (1971) Appl. Microbiol . 22. 17-22.
11. Chadwick P., Delisle G.H. and Byer M. (1974) Can. J. Microbiol . 20. 1653-1664.

Composition:

SS AGAR:

Lactose - 10.0 g/L
Meat Peptone - 2.5 g/L
Bile Salts - 8.5 g/L
Ferric Citrate - 1.0 g/L
Sodium Citrate - 8.5 g/L
Neutral Red - 25.0 mg/L
Sodium Thiosulfate - 8.5 g/L
Brilliant Green - 0.33 mg/L
Beef Extract - 5.0 g/L
Agar - 13.5 g/L
Casein Peptone - 2.5 g/L

XLD AGAR:

Yeast extract - 3.0 g/L
L-Lysine HCl - 5.0 g/L
Xylose - 3.75 g/L
Lactose - 7.5 g/L
Sucrose - 7.5 g/L
Sodium desoxycholate - 1.0 g/L
Sodium chloride - 5.0 g/L
Sodium thiosulphate - 6.8 g/L
Ferric ammonium citrate - 0.8 g/L
Phenol red - 0.08 g/L
Agar - 12.5 g/L

Storage: 2-8 °C

Appearance: SS AGAR: Red-Orange, Very Slightly Opalescent

XLD AGAR: Red coloured

pH Range: SS AGAR: 6.8 - 7.2

XLD AGAR: 7.2 – 7.6

Package contents: 10 plates in a package

Exp. Date: Printed on label and on the item.

Required materials not supplied: Laboratory equipment as required.
Hy Laboratories Ltd.

6 Menachem Plaut St., Park Tamar, Rehovot 7670606, Israel
Tel. +972.8.9366475 Email. hylabs@hylabs.co.il

Implementation Date: 28/10/25
Version Number: 01

Warning and Precautions - For professional use only. Follow good microbiological lab practices while handling specimens and culture. Do not use Petri dishes if they show evidence of microbial contamination, discoloration, drying, cracking, or other signs of deterioration. Avoid freezing and overheating. The Petri Dishes may be used / inoculated up to the expiration date and incubated for the recommended incubation times. After use and prior to discarding, specimen containers and all contaminated material, including the used culture media and contaminated culture containers, must be sterilized or incinerated by validated procedures. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.

If excessive moisture is observed, invert the bottom over an off-set lid and allow to air dry in order to prevent formation of a seal between the top and bottom of the plate during incubation. Storage Instructions: On receipt, store plates in the dark at 2–8°C. Avoid freezing and overheating. Do not open until ready to use.

Performance Testing Results:

GPT: inoculum 10-100 cfu.

Inhibitory properties: inoculum 10000 cfu.

TEST	ATCC	Incubation Temp. (°C)	Incubation Cond.	Reaction 1 SS AGAR		Reaction 2 XLD AGAR	
<i>Salmonella typhimurium</i>	14028	33-37 °C	Aerobic, 18-24 hours	Growth	Colorless, black centered	Growth	Pink-red, black centered
<i>Shigella flexneri</i>	29903	33-37 °C	Aerobic, 18-24 hours	Growth	Smooth opaque colorless	Growth	Pink-red
<i>Shigella sonnei</i>	29930	33-37 °C	Aerobic, 18-24 hours	Growth	Colorless	Growth	Smooth, opaque, colorless
<i>Proteus mirabilis</i>	4630	33-37 °C	Aerobic, 18-24 hours	Growth	Smooth, opaque, may be black	Growth	Pink to orange, may be black
<i>Escherichia coli</i>	25922	33-37 °C	Aerobic, 18-24 hours	Poor	Pink	Poor	Yellow
<i>Enterococcus faecalis</i>	19433	33-37 °C	Aerobic, 18-24 hours	Inhibition		Inhibition	