

Bioneer's AccuPower™ PCR PreMix

*When accuracy is the question....
AccuPower is the answer*

Add Template & Primers



K- 2012	AccuPower PCR PreMix. Only 139NIS for 96 PreMix reactions.	Buy 5 for 556NIS (K-2016) 1.15NIS per reaction
K- 5050	AccuPower HotStart PCR PreMix. Only 169NIS for 96 PreMix reactions.	Buy 5 for 676NIS (K-5051) 1.4NIS per reaction
K- 3201	Accupower TLA PCR PreMix. High fidelity, High yield with blunt end. Only 235NIS for 96 PreMix reactions.	Buy 5 for 940NIS (K-3203) 1.95NIS per reaction
K- 3301	Accupower HF PCR PreMix. High fidelity (HotStart), High yield with A end Only 245NIS for 96 PreMix reactions.	Buy 5 for 980NIS (K-3303) 2NIS per reaction
K- 2041	AccuPower RT PreMix. Only 396NIS for 96 PreMix reactions.	
K- 2055	AccuPower RT/PCR PreMix Only 576NIS for 96 PreMix reactions (one step).	

- Accupower PreMix contain 96 ready to use tubes (just add template, primers and DDW)
- The prices are for 0.2ml tubes and for 20µl rxn volume.
- Accupower is stable at -20°C for a few years.

For more information and a free sample please contact:

Shachar: 052-8343519 shachar@hylabs.co.il