



Hy-Laboratories Ltd.
Park Tamar
Rehovot 76326 Israel
Tel : 972-8-9366475
Fax : 972-8-9366474
hylabs@hylabs.co.il
www.hylabs.co.il

YERSINIAE SELECTIVE AGAR

Instructions for use

Cat No. PD097

Inoculate with approximately 10^2 c.f.u. of the test cultures and incubate at 22-32°C. Examine the plates for growth after 18-24 hours and 48 hours.

RESULTS

Organisms	Growth	Colony Morphology
<i>Yersinia enterocolitica</i>	Good	Translucent with dark pink center*
<i>E. coli</i>	Inhibited	
<i>Proteus mirabilis</i>	Inhibited	
<i>Pseudomonas aeruginosa</i>	Inhibited	
<i>Strep. faecalis</i>	Inhibited	

* Bile precipitation may be present on the agar.