



Hy-Laboratories Ltd.  
Park Tamar  
Rehovot 76326 Israel  
Tel : 972-8-9366475  
Fax : 972-8-9366474  
[hylabs@hylabs.co.il](mailto:hylabs@hylabs.co.il)

## CDC Anaerobic Blood Agar

### Cat. No. PD162

Culture media is a crucial factor in the success of the anaerobic culture. CDC Agar is supplemented with sheep blood and is presently recognized as the best non-selective solid medium for growth of anaerobes due to the numerous growth factors including: L-cystine, Hemin and Vit K1. It is recommended for primary isolation of all types of anaerobes found in clinica, including those requiring additional growth factors, such as *Bac. melaninogenicus*, *F. necrophorum*, etc. which grow with typical morphological characteristic colonies on this medium.