

Fluid A

A rinsing solution

BP316 Fluid A, different volume bottle with screw cap or septum cap

Appearance: Clear

Intended Use and Principles of the Procedure – Fluid A is a rinsing solution, suitable for sterility testing or for bioburden testing using membrane filtration in compliance with the United States Pharmacopeia (USP) and European Pharmacopeia (EP).

Fluid A contains 0.1% peptone. Peptone is used as a nutritive substrate for different microorganisms.

Formula

Peptone No. 3 - 1 gr/Liter

QC performance

pH (adjusted) at RT - 7.1 ± 0.2

Storage 15-25 °C

Warning and Precautions - For laboratory use only. Do not use bottles if they show evidence of microbial contamination, discoloration, cracking, or other signs of deterioration. Do not use the product beyond its expiry date