



Hy-Laboratories Ltd.  
Park Tamar  
Rehovot 76326 Israel  
Tel : 972-8-9366475  
Fax : 972-8-9366474  
[hylabs@hylabs.co.il](mailto:hylabs@hylabs.co.il)  
[www.hylabs.co.il](http://www.hylabs.co.il)

## **BOTTLED PRODUCTS (BP\*) – for food testing**

---

<b>Cat. No</b>	<b>Type of Media</b>	<b>Recommended use</b>
BP202	APT AGAR	Isolation and enumeration of heterofermentative lactobacilli and lactic streptococci.
BP297	BAIRED PARKER AGAR BASE	Baired parker agar base is used with egg yolk tellurite** or RPF** for selective isolation and enumeration of coagulase- positive staphlococci .
BP262	CHROMagar E.Coli (improved TBX)	Isolation and differentiation of E. coli by specific chromogenic enzymatic reaction.
BP207	OGY AGAR	Enumeration of Yeasts and Molds. Note: Add ox tetracycline** Immediately before pouring for Selective, inhibits growth of bacteria
BP442	ORANGE SERUM AGAR	Cultivation and enumeration of microorganism associated with the spoilage of citruss products.
BP032	SMA	Total count of microorganisms as recommended by HPHA.
BP210	VRBA	Selective medium for detection, enumeration and identification of coliform.
BP358	VRBG	Detection and enumeration of enterobacteriaceae.

---

\*avaibale in different volum.

\*\*Also available upon request.