



Hy-Laboratories Ltd.  
Park Tamar  
Rehovot 76326 Israel  
Tel : 972-8-9366475  
Fax : 972-8-9366474  
[hylabs@hylabs.co.il](mailto:hylabs@hylabs.co.il)  
[www.hylabs.co.il](http://www.hylabs.co.il)

## CHROMagar S. Aureus

**PD235 , DD066 , DS049 , DS050**

Chromogenic medium for the isolation of *Staphylococcus aureus*.

A time saving breakthrough in microbiological testing:  
**Streak, Inoculate, Identify.**

A color detection of *S. aureus* based on it's enzymatic activity on chromogenic substracts.

*S. aureus* gives mauve-pink colonies after 24 hours. Other microorganisms are inhibited or grow in blue, colorless or cream clones.

The divided plate DD044, CHROMagar S. aureus and CHROMagar MRSA, with antibiotics addition at a final concentration of 4 mg cefoxitim litse/L, allows the isolation of methicillin resistant staphylococci.