



Hy-Laboratories Ltd.
Park Tamar
Rehovot 76326 Israel
Tel : 972-8-9366475
Fax : 972-8-9366474
hylabs@hylabs.co.il
www.hylabs.co.il

CHROMagar E. Coli

LD503, BP262, DS061, DS059

Chromogenic and selective media for the detection ,
identification and enumeration of *Escherichia coli* .

Highly selective * Bile salts : (Inhibition of Gram positive
Bacteria)

Differential * X- GLUC substrate : (5-bromo-4 chloro-
3 indolyl β - D glucuronide.)

β – glucuronidase positive (green blue colonies): *E. coli*

β – glucuronidase negative (colorless colonies) : Other bacteria