

CARY-BLAIR (CB) TRANSPORT MEDIUM

Cat No. TT225

Proper collection and preservation of feces is an important requirement

for the isolation of microorganisms producing the following intestinal disease: *Shigella*, *Salmonella*.

Campylobacter, *Vibrio sp.* When unavoidable delays occur in transporting the specimen to the laboratory, these microorganisms will not survive as the changes in pH occurs with a drop in temperature.

The Cary-Blair medium maintains the microorganism's viability during transport and supports a better recovery of *Shigella* than Stuart's medium, and recovery of *Cholera*, *Salmonella* and *Yersinia* after several days. This medium is available at Hy-Labs in the following formulation:

Plastic vial and Spork containing Cary Blair (CB) Medium.

For transport of *Shigella*, *Salmonella*, *Yersinia*, *Vibrio spp.*