



Hy-Laboratories Ltd.
Park Tamar
Rehovot 76326 Israel
Tel : 972-8-9366475
Fax : 972-8-9366474
hylabs@hylabs.co.il
www.hylabs.co.il

CHROMagar S. Aureus

PD235 , DD066 , DS049 , DS050

Chromogenic medium for the isolation of *Staphylococcus aureus*.

A time saving breakthrough in microbiological testing:
Streak, Inoculate, Identify.

A color detection of *S. aureus* based on it's enzymatic activity on chromogenic substracts.

S. aureus gives mauve-pink colonies after 24 hours. Other microorganisms are inhibited or grow in blue, colorless or cream colonies.

The divided plate DD044, CHROMagar S. aureus and CHROMagar MRSA, with antibiotics addition at a final concentration of 4 mg cefoxitim litse/L, allows the isolation of methicillin resistant staphylococci.