


CHROMagar™ LIM QUICK Strep B Broth

TT – 390

For rapid detection of *Streptococcus agalactiae* from swabs.

- LIM RambaQUICK™ StrepB is a selective enrichment broth designed to largely reduce the growth of Enterococci, while allowing growth of GBS, and, therefore to raise the sensitivity of the prepartum screening method.

LIM QUICK™ StrepB Method



Performance of the LIM QUICK™ StrepB broth

- Compared to classical LIM + subculture method:
INHIBITION OF MOST ENTEROCOCCI
LIM QUICK™ StrepB is a selective enrichment broth designed, and, therefore to raise the sensitivity the prepartum screening method.
- **BROTH INOCULATION**
Classical Method recommends overnight broth incubation for classical LIM.
LIM RambaQUICK™ StrepB can be used also with overnight incubation, and has also an optimal differentiation at 6h only.